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HiScribe™ T7 mRNA 合成试剂盒 (含 CleanCap® Reagent AG)

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- Dust & air
- Skin & hair
- Aqueous solutions & reagents
- Most surfaces (doorknobs, keyboards)

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- Wear laboratory gloves and change them often
- Use RNase-free certified, disposable plasticware and solutions
- Decontaminate glassware & plasticware
- Maintain a separate, clean surface for RNA work

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Diethylpyrocarbonate (DEPC) treatment:

1. Add 1 ml DEPC per liter of solution.
2. Stir for 1 hour.
3. Autoclave for 1 hour.

⊗ DEPC:

- Compounds with primary amine groups (e.g., Tris)
- Compounds that are not stable during autoclaving

Not using DEPC?

Prepare solution with Nuclease-free Water or Milli-Q® water

Dissolving solids:

- Use high-purity solids (e.g., DTT, nucleotides, manganese salts)
- Use autoclaved DEPC-treated or Milli-Q water
- Sterilize with a 0.22 µm filter

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模板生成	体外转录		RNA 加帽	Poly (A) 加尾	RNA 纯化
Q5® High-Fidelity DNA Polymerase	HiScribe® T7 mRNA Kit with CleanCap® Reagent AG			Monarch® RNA Cleanup Kit (10 µg)	
	HiScribe T7 ARCA mRNA Synthesis Kit (with tailing)				
dNTP solution mixes	HiScribe T7 ARCA mRNA Synthesis Kit		E. coli Poly(A) Polymerase	Monarch RNA Cleanup Kit (50 µg)	
GMP COMING SOON BspQI*	GMP HiScribe T7 High Yield RNA Synthesis Components	GMP Vaccinia Capping System		Monarch RNA Cleanup Kit (500 µg)	
DNA Assembly - NEBuilder HiFi DNA Assembly - Golden Gate Assembly	HiScribe T7 Quick High Yield RNA Synthesis Kit	GMP mRNA Cap 2'-O-Methyltransferase		Lithium Chloride	
	HiScribe SP6 High Yield RNA Synthesis Kit	ARCA and other mRNA cap analogs			
	T3 & SP6 RNA Polymerases				
	GMP T7 RNA Polymerase				
	Hi-T7 RNA Polymerase				
	Companion Products				
	GMP RNase inhibitor (Murine)			Monarch RNA Cleanup Binding Buffer	
	RNase Inhibitor (Human Placental)			Monarch RNA Cleanup Wash Buffer	
	GMP Pyrophosphatase, Inorganic (E. coli)			Nuclease-free Water	
	Pyrophosphatase, Inorganic (Yeast)				
	GMP DNase I (RNase-free)				
	NTPs				

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产品名称	货号	规格
HiScribe T7 mRNA Kit with CleanCap Reagent AG	E2080S	20 reactions
HiScribe T7 Quick High Yield RNA Synthesis Kit	E2050S	50 reactions
HiScribe T7 High Yield RNA Synthesis Kit	E2040S	50 reactions
HiScribe SP6 RNA Synthesis Kit	E2070S	50 reactions
HiScribe T7 ARCA mRNA Kit	E2065S	20 reactions
HiScribe T7 ARCA mRNA Kit (with tailing)	E2060S	20 reactions
T3 RNA Polymerase	M0378S	5,000 units
T7 RNA Polymerase	M0251S/L	5,000/25,000 units
SP6 RNA Polymerase	M0207S/L	2,000/10,000 units
Hi-T7 RNA Polymerase	M0658S	5,000 units
<i>E. coli</i> Poly(A) Polymerase	M0276S/L	100/500 units
Poly(U) Polymerase	M0337S	60 units
<i>E. coli</i> RNA Polymerase, Core Enzyme	M0550S	100 units
<i>E. coli</i> RNA Polymerase, Holoenzyme	M0551S	50 units
Ribonucleotide Solution Set	N0450S/L	10/50 µmol of each
Ribonucleotide Solution Mix	N0466S/L	10/50 µmol of each
Pyrophosphatase, Inorganic (<i>E. coli</i>)	M0361S/L	10/50 units
Pyrophosphatase, Inorganic (yeast)	M2403S/L	10/50 units
Thermostable Inorganic Pyrophosphatase	M0296S/L	250/1,250 units
mRNA Decapping Enzyme	M0608S	2,000 units
Vaccinia Capping System	M2080S	400 units
Anti-Reverse Cap Analog 3'-O-Me-m ⁷ G(5')ppp(5')G	S1411S/L	1/5 µmol
Standard Cap Analog m ⁷ G(5')ppp(5')G	S1404S/L	1/5 µmol
Unmethylated Cap Analog G (5')ppp(5')G	S1407S/L	1/5 µmol
Methylated Cap Analog for A + 1 sites m ⁷ G(5')ppp(5')A	S1405S/L	1/5 µmol
Unmethylated Cap Analog for A + 1 sites G(5')ppp(5')A	S1406S/L	1/5 µmol
mRNA Cap 2'-O-Methyltransferase	M0366S	2,000 units
3'-Desthiobiotin-GTP	N0761S	0.5 µmol

cDNA 合成

产品名称	货号	规格
ProtoScript II Reverse Transcriptase	M0368S/L/X	4,000/10,000/40,000 units
M-MuLV Reverse Transcriptase	M0253S/L	10,000/50,000 units
AMV Reverse Transcriptase	M0277S/L	200/1,000 units
WarmStart RTx Reverse Transcriptase	M0380S/L	50/250 reactions
ProtoScript II First Strand cDNA Synthesis Kit	E6560S/L	30/150 reactions
ProtoScript First Strand cDNA Synthesis Kit	E6300S/L	30/150 reactions

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产品名称	货号	规格
Luna Universal qPCR Master Mix	M3003S/L/X/E	200/500/1,000/2,500 reactions
Luna Universal Probe qPCR Master Mix	M3004S/L/X/E	200/500/1,000/2,500 reactions
Luna Universal One-Step RT-qPCR Kit	E3005S/L/X/E	200/500/1,000/2,500 reactions
Luna Universal Probe One-Step RT-qPCR Kit	E3006S/L/X/E	200/500/1,000/2,500 reactions

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产品名称	货号	规格
Luna Probe One-Step RT-qPCR Kit (No ROX)	E3007E	2,500 reactions
LunaScript RT SuperMix Kit	E3010S/L	25/100 reactions
LunaScript RT Master Mix Kit (Primer-free)	E3025S/L	25/100 reactions
Luna Probe One-Step RT-qPCR 4X Mix with UDG	M3019S/L/X/E	200/500/1,000/2,500 reactions
Luna SARS-CoV-2 RT-qPCR Multiplex Assay Kit	E3019S/L	96/480 reactions

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产品名称	货号	规格
T4 RNA Ligase 1 (ssRNA Ligase)	M0204S/L	1,000/5,000 units
T4 RNA Ligase 1 (ssRNA Ligase), high concentration	M0437M	5,000 units
T4 RNA Ligase 2 (dsRNA Ligase)	M0239S/L	150/750 units
T4 RNA Ligase 2, truncated	M0242S/L	2,000/10,000 units
T4 RNA Ligase 2, truncated K227Q	M0351S/L	2,000/10,000 units
T4 RNA Ligase 2, truncated KQ	M0373S/L	2,000/10,000 units
RtcB Ligase	M0458S/L	25 reactions
Thermostable 5' AppDNA/RNA Ligase	M0319S/L	10/50 reactions
5' DNA Adenylation Kit	E2610S/L	10/50 reactions
SplintR Ligase	M0375S/L	1,250/6,250 units
RNA 5' Pyrophosphohydrolase (RppH)	M0356S	200 units
5' Deadenylase	M0331S	1,000 units
RNase I ₁	M0243S/L	5,000/25,000 units
RNase H	M0297S/L	250/1,250 units
RNase HII	M0288S/L	250/1,250 units
Quick Dephosphorylation Kit	M0508S/L	100/500 units
Antarctic Phosphatase	M0289S/L	1,000/5,000 units
Alkaline Phosphatase Calf Intestinal (CIP)	M0290S/L	1,000/5,000 units
Shrimp Alkaline Phosphatase (rSAP)	M0371S/L	500/2,500 units
T4 Polynucleotide Kinase	M0201S/L	500/2,500 units
ShortCut RNase III	M0245S/L	200/1,000 units
XRN-1	M0338S/L	20/100 units
Exonuclease T	M0265S/L	250/1,250 units
NudC Pyrophosphatase	M0607	250 pmol
Sce PUS1	M0526	5,000 pmol

gRNA 合成

产品名称	货号	规格
EnGen sgRNA Synthesis Kit, <i>S. pyogenes</i>	E3322S	20 reactions

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产品名称	货号	规格
RNase Inhibitor, Murine	M0314S/L	3,000/15,000 units
RNase Inhibitor, Human Placenta	M0307S/L	2,000/10,000 units
Ribonucleoside Vanadyl Complex	S1402S	10 ml (200 mM)

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产品名称	货号	规格
NEBNext Single Cell/Low Input RNA Library Prep Kit for Illumina	E6420S/L	24/96 reactions
NEBNext Single Cell/Low Input cDNA Synthesis & Amplification Module	E6421S/L	24/96 reactions
NEBNext Ultra II Directional RNA Library Prep Kit for Illumina	E7760S/L	24/96 reactions
NEBNext Ultra II Directional RNA Library Prep with Sample Purification Beads	E7765S/L	24/96 reactions
NEBNext Ultra II RNA Library Prep Kit for Illumina	E7770S/L	24/96 reactions
NEBNext Ultra II RNA Library Prep with Sample Purification Beads	E7775S/L	24/96 reactions
NEBNext Multiplex Small RNA Library Prep Set for Illumina (Set 1)	E7300S/L	24/96 reactions
NEBNext Multiplex Small RNA Library Prep Set for Illumina (Set 2)	E7580S/L	24/96 reactions
NEBNext Multiplex Small RNA Library Prep Kit for Illumina (Index Primers 1-48)	E7560S	96 reactions
NEBNext Small RNA Library Prep Set for Illumina (Multiplex Compatible)	E7330S/L	24/96 reactions
NEBNext rRNA Depletion Kit v2 (Human/Mouse/Rat)	E7400S/L/X	6/24/96 reactions
NEBNext rRNA Depletion Kit v2 (Human/Mouse/Rat) with RNA Sample Purification Beads	E7405S/L/X	6/24/96 reactions
NEBNext rRNA Depletion Kit (Bacteria)	E7850S/L/X	6/24/96 reactions
NEBNext rRNA Depletion Kit (Bacteria) with RNA Sample Purification Beads	E7860S/L/X	6/24/96 reactions
NEBNext Globin & rRNA Depletion Kit (Human/Mouse/Rat)	E7750S/L/X	6/24/96 reactions
NEBNext Globin & rRNA Depletion Kit (Human/Mouse/Rat) with RNA Sample Purification Beads	E7755S/L/X	6/24/96 reactions
NEBNext RNA Depletion Core Reagent Set	E7865S/L/X	6/24/96 reactions
NEBNext RNA Depletion Core Reagent Set with RNA Sample Purification Beads	E7870S/L/X	6/24/96 reactions
NEBNext Poly(A) mRNA Magnetic Isolation Module	E7490S/L	24/96 reactions
NEBNext Library Quant Kit for Illumina	E7630S/L	100/500 reactions
NEBNext Multiplex Oligos for Illumina (Index Primers Set 1)	E7335S/L	24/96 reactions
NEBNext Multiplex Oligos for Illumina (Index Primers Set 2)	E7500S/L	24/96 reactions
NEBNext Multiplex Oligos for Illumina (Index Primers Set 3)	E7710S/L	24/96 reactions
NEBNext Multiplex Oligos for Illumina (Index Primers Set 4)	E7730S/L	24/96 reactions
NEBNext Multiplex Oligos for Illumina (96 Index Primers)	E6609S/L	96/384 reactions
NEBNext Multiplex Oligos for Illumina (Dual Index Primers Set 1)	E7600S	96 reactions

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产品名称	货号	规格
Monarch Total RNA Miniprep Kit	T2010S	50 preps
Monarch RNA Lysis Buffer	T2012L	100 ml
Monarch RNA Priming Buffer	T2013L	56 ml
Monarch RNA Wash Buffer	T2014L	50 ml
Monarch Collection Tubes II	T2018L	100 tubes
Monarch RNA Purification Columns	T2007L	100 columns
Monarch DNA/RNA Protection Reagent	T2011L	56 ml
Monarch Total RNA Miniprep Enzyme Pack (contains DNase I, Prot K, and associated buffers)	T2019L	1 pack
Monarch RNA Cleanup Kit (10 µg)	T2030S/L	10/100 preps
Monarch RNA Cleanup Kit (50 µg)	T2040S/L	10/100 preps
Monarch RNA Cleanup Kit (500 µg)	T2050S/L	10/100 preps
Monarch RNA Cleanup Binding Buffer	T2041L	80 ml
Monarch RNA Cleanup Wash Buffer	T2042L	40 ml
Monarch RNA Cleanup Columns (10 µg)	T2037L	100 columns and tubes
Monarch RNA Cleanup Columns (50 µg)	T2047L	100 columns and tubes
Monarch RNA Cleanup Columns (500 µg)	T2057L	100 columns and tubes
Magnetic mRNA Isolation Kit	S1550S	25 isolations
Epimark N6-Methyladenosine Enrichment Kit	E1610S	20 reactions
Oligo d(T) ₂₅ Magnetic Beads	S1419S	25 mg
Oligo d(T) ₂₅ Cellulose Beads	S1408S	250 mg
Streptavidin Magnetic Beads	S1420S	5 ml (20 mg)
Hydrophilic Streptavidin Magnetic Beads	S1421S	5 ml (20 mg)
polyA Spin mRNA Isolation Kit	S1560S	8 isolations
p19 siRNA Binding Protein	M0310S	1,000 units

RNA Markers & Ladders

产品名称	货号	规格
dsRNA Ladder	N0363S	25 gel lanes
microRNA Marker	N2102S	100 gel lanes
ssRNA Ladder	N0362S	25 gel lanes
Low Range ssRNA Ladder	N0364S	25 gel lanes
RNA Loading Dye (2X)	B0363S	4 ml
Universal miRNA Cloning Linker	S1315S	0.83 nmol

NEB (北京)
客服及技术支持
客服电话: 400-811-2220/400-690-3366
电话: 010-82378266
传真: 010-82378262
e-mail: support@neb-china.com
网址: www.neb-china.com
官方微信: NEBiolabs

美国
New England Biolabs, Inc.
Telephone (978) 927-5054
Toll Free (USA Orders) 1-800-632-5227
Toll Free (USA Tech) 1-800-632-7799
info@neb.com

澳大利亚 / 新西兰
New England Biolabs (Australia) PTY
Telephone: 1800 934 218 (AU)
info.au@neb.com
Telephone: 0800 437 209 (NZ)
info.nz@neb.com

加拿大
New England Biolabs, Ltd.
Toll Free: 1-800-387-1095
info.ca@neb.com

法国
New England Biolabs France
Telephone: 0800 100 632
info.fr@neb.com

德国 / 奥地利
New England Biolabs GmbH
Free Call: 0800/246 5227
(Germany)
Free Call: 00800/246 52277
(Austria)
info.de@neb.com

日本
New England Biolabs Japan, Inc.
Telephone: +81 (0)3 5669 6191
info.jp@neb.com

新加坡
New England Biolabs, Pte. Ltd.
Telephone: +65 638 59623
sales.sg@neb.com

英国
New England Biolabs (UK), Ltd.
Call Free: 0800 318486
info.uk@neb.com

www.neb.com

One or more of these products are covered by patents, trademarks and/or copyrights owned or controlled by New England Biolabs, Inc. For more information, please email us at busdev@neb.com. The use of these products may require you to obtain additional third party intellectual property rights for certain applications.

Your purchase, acceptance, and/or payment of and for NEB's products is pursuant to NEB's Terms of Sale at www.neb.com/support/terms-of-sale. NEB does not agree to and is not bound by any other terms or conditions, unless those terms and conditions have been expressly agreed to in writing by a duly authorized officer of NEB.

B CORPORATION™ is a trademark of B Lab, Inc. CLEANCAP® is a registered trademark of TriLink BioTechnologies, LLC. CLONTECH® is a registered trademark of Clontech Laboratories. SMART-SEQ® is a registered trademark of Illumina, Inc. NEXTERA® is a registered trademark of FPL Group, Inc. KAPA® is a registered trademark of Kapa Biosystems. ZYMOLYASE® is a registered trademark of Zymo Research. CLEAN AND CONCENTRATOR™ is a trademark of Zymo Research. SYBR® is a registered trademark of Molecular Probes, Inc. AGILENT® and BIOANALYZER® is a registered trademark of Agilent Technologies, Inc. TAQMAN® and NANODROP® are a registered trademarks of ThermoFisher Scientific, Inc. TRIZOL® is a registered trademark of Molecular Research Center, Inc. ILLUMINA®, MISEQ® and NEXTSEQ® a registered trademark of Illumina, Inc.

© Copyright 2021, New England Biolabs, Inc.; all rights reserved

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RNA_TG – Version 4.0 – 10/21

Download the **NEB AR App** for iPhone® or iPad®. Scan the augmented reality butterfly icon located on the corner of the page to find videos, tutorials and immersive experiences.

RNA Markers & Ladders