
New England Biolabs Product Specification

Product Name: Rapid PNGase F Antibody Standard
Catalog #: P6043S
Concentration: 5 mg/ml
Unit Definition: N/A
Shelf Life: 24 months
Storage Temp: -20°C
Storage Conditions: 100 mM NaCl, 10 mM Sodium Phosphate, (pH 7.0 @ 25°C)
Specification Version: PS-P6043S v1.0
Effective Date: 23 Nov 2015

Assay Name/Specification (minimum release criteria)

Functional Testing (Rapid Deglycosylation) - A 10 µl reaction in 1X Rapid PNGase F Buffer containing 32 µg of Rapid PNGase F Antibody Standard and 1 µl Rapid PNGase F incubated for 10 minutes at 50°C results in a mobility shift in > 95% of the starting material as determined by Coomassie stained SDS-PAGE.



Date 23 Nov 2015

Derek Robinson
Director of Quality Control

