

---

## New England Biolabs Product Specification

*Product Name:* Anti-SNAP-tag<sup>®</sup> Antibody (Polyclonal)  
*Catalog #:* P9310S  
*Shelf Life:* 24 months  
*Storage Temp:* -20°C  
*Storage Conditions:* 150 mM NaCl, 10 mM HEPES, 50 % Glycerol, (pH 7.5 @ 25°C)  
*Specification Version:* PS-P9310S v1.0  
*Effective Date:* 05 Feb 2016

<b>Assay Name/Specification (minimum release criteria)</b>
------------------------------------------------------------

<b>Functional Testing (Western Blot, Antibody)</b> - A 1:1,000 dilution of the antibody is tested for its ability to detect 5 ng of SNAP-tag purified Protein by Western Blot.
--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------



Date 05 Feb 2016

---

Derek Robinson  
Director of Quality Control

