
New England Biolabs Product Specification

Product Name: Monarch[®] Gel Dissolving Buffer
Catalog #: T1021L
Shelf Life: 24 months
Storage Temp: 25°C
Composition: Proprietary
Specification Version: PS-T1021L v1.0
Effective Date: 08 Nov 2017

Assay Name/Specification (minimum release criteria)

Functional Testing (Gel Purification and Ligation) - A linear pLA3 fragment is gel purified and added to 8 Monarch [®] DNA Cleanup Columns, eluted and ligated with Blunt/TA Ligase Master Mix resulting in ≥90% conversion to a multiband pattern observed by gel electrophoresis in ≥75% of the samples.
--



Date 08 Nov 2017

Derek Robinson
Director of Quality Control

